

For Research Use Only

FAM168B Polyclonal antibody

Catalog Number: 24422-1-AP



Basic Information

Catalog Number:

24422-1-AP

Size:

850 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19535

GenBank Accession Number:

BC066347

GeneID (NCBI):

130074

UNIPROT ID:

A1KXE4

Full Name:

family with sequence similarity 168, member B

Calculated MW:

195 aa, 20 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:20-1:200

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC : human cervical cancer tissue, human stomach cancer tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human cervical cancer slide using 24422-1-AP (FAM168B Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human cervical cancer slide using 24422-1-AP (FAM168B Antibody) at dilution of 1:50 (under 40x lens).