

For Research Use Only

# FAM168B Polyclonal antibody

Catalog Number: 24422-1-AP



## Basic Information

<b>Catalog Number:</b> 24422-1-AP	<b>GenBank Accession Number:</b> BC066347	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 850 µg/ml	<b>GeneID (NCBI):</b> 130074	<b>Recommended Dilutions:</b> IHC 1:20-1:200
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> A1KXE4	
<b>Isotype:</b> IgG	<b>Full Name:</b> family with sequence similarity 168, member B	
<b>Immunogen Catalog Number:</b> AG19535	<b>Calculated MW:</b> 195 aa, 20 kDa	

## Applications

**Tested Applications:**  
IHC, ELISA

**Species Specificity:**  
human

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

**Positive Controls:**

IHC: human cervical cancer tissue, human stomach cancer tissue

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human cervical cancer slide using 24422-1-AP (FAM168B Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human cervical cancer slide using 24422-1-AP (FAM168B Antibody) at dilution of 1:50 (under 40x lens).