

For Research Use Only

LRP4 Polyclonal antibody, PBS Only

Catalog Number: 24434-1-PBS



Basic Information

Catalog Number:

24434-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19487

GenBank Accession Number:

BC136667

GeneID (NCBI):

4038

UNIPROT ID:

O75096

Full Name:

low density lipoprotein receptor-related protein 4

Calculated MW:

1905 aa, 212 kDa

Observed MW:

240 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

LRP4 (low-density lipoprotein receptor-related protein 4), also known as multiple epidermal growth factor-like domains 7 (MEGF7), is a type I single transmembrane protein of the low-density lipoprotein receptor (LDLR) family. LRP4 is composed of a large extracellular region containing multiple LDLa repeats, EGF-like domains and β -propeller domains, a transmembrane region and a cytoplasmic region containing an NPxY motif and a PDZ-interacting motif (PMID: 26071838). LRP4 binds agrin and recruits it to the MuSK signaling complex, which is critical for neuromuscular junction formation (PMID: 17119023; 18957220). LRP4 mutation or malfunction is implicated in disorders including congenital myasthenic syndrome, myasthenia gravis, and diseases of bone or kidney (26071838).

Storage

Storage:

Store at -80°C .

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Mouse brain tissue was subjected to SDS PAGE followed by western blot with 24434-1-AP (LRP4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24434-1-PBS in a different storage buffer formulation.