

For Research Use Only

FBXL17 Polyclonal antibody

Catalog Number: 24476-1-AP



Basic Information

Catalog Number:

24476-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20436

GenBank Accession Number:

BC126144

GeneID (NCBI):

64839

UNIPROT ID:

Q9UF56

Full Name:

F-box and leucine-rich repeat protein
17

Calculated MW:

701 aa, 76 kDa

Observed MW:

76 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB: mouse testis tissue, rat testis tissue

Background Information

FBXL17 (F-box and leucine-rich repeat protein 17) is a substrate-recognition component of the SCF-type E3 ubiquitin ligase complex (SKP1-CUL1-F-box).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24476-1-AP (FBXL17 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.