

For Research Use Only

TMEM194A Polyclonal antibody

Catalog Number: 24477-1-AP



Basic Information

Catalog Number:

24477-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21614

GenBank Accession Number:

BC126300

GeneID (NCBI):

23306

UNIPROT ID:

O14524

Full Name:

transmembrane protein 194A

Calculated MW:

444 aa, 51 kDa

Observed MW:

40-42 kDa, 50-55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : MCF-7 cells, mouse liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF-7 cells were subjected to SDS PAGE followed by western blot with 24477-1-AP (TMEM194A Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.