

For Research Use Only

DNAH17 Polyclonal antibody, PBS Only

Catalog Number: 24488-1-PBS



Basic Information

Catalog Number:

24488-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20033

GenBank Accession Number:

BC112284

GeneID (NCBI):

8632

UNIPROT ID:

Q9UFH2

Full Name:

dynein, axonemal, heavy chain 17

Calculated MW:

4485 aa, 512 kDa

Observed MW:

500-512 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to Tris-Acetate gel system followed by western blot with 24488-1-AP (DNAH17 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24488-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 24488-1-AP (DNAH17 Antibody) at dilution of 1:50 (under 10x lens). This data was developed using the same antibody clone with 24488-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 24488-1-AP (DNAH17 Antibody) at dilution of 1:50 (under 40x lens). This data was developed using the same antibody clone with 24488-1-PBS in a different storage buffer formulation.