

For Research Use Only

LCK Polyclonal antibody

Catalog Number: 24497-1-AP

1 Publications



Basic Information

Catalog Number:

24497-1-AP

Concentration:

950 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21475

GenBank Accession Number:

BC013200

GeneID (NCBI):

3932

UNIPROT ID:

P06239

Full Name:

Lymphocyte-specific protein tyrosine kinase

Calculated MW:

539 aa, 56 kDa

Observed MW:

50-60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse

Positive Controls:

WB: mouse thymus tissue, Jurkat cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 24497-1-AP (LCK Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.