

For Research Use Only

PPP2R3B Polyclonal antibody, PBS Only

Catalog Number: 24502-1-PBS



Basic Information

Catalog Number:

24502-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21366

GenBank Accession Number:

BC009032

GeneID (NCBI):

28227

UNIPROT ID:

Q9Y5P8

Full Name:

protein phosphatase 2 (formerly 2A),
regulatory subunit B^β, beta

Calculated MW:

65 kDa

Observed MW:

65 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

PPP2R3B stands for Protein Phosphatase 2 Regulatory Subunit B^β beta. It is a protein-coding gene that plays a vital role in regulating cellular processes. Gonosome-linked expression of PPP2R3B in cutaneous melanoma correlates with distant metastasis free survival. It encodes the PR70 protein (also known as PR48), which acts as a regulatory subunit for Protein Phosphatase 2A (PP2A)

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HL-60 cells were subjected to SDS PAGE followed by western blot with 24502-1-AP (PPP2R3B Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24502-1-PBS in a different storage buffer formulation.