

For Research Use Only

C6orf173 Polyclonal antibody

Catalog Number: 24509-1-AP



Basic Information

Catalog Number:

24509-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21580

GenBank Accession Number:

BC046178

GeneID (NCBI):

387103

UNIPROT ID:

Q5EE01

Full Name:

chromosome 6 open reading frame

173

Calculated MW:

88 aa, 10 kDa

Observed MW:

7-10 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : L02 cells, A549 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



L02 cells were subjected to SDS PAGE followed by western blot with 24509-1-AP (C6orf173 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.