

For Research Use Only

TBC1D13 Polyclonal antibody

Catalog Number: 24562-1-AP



Basic Information

Catalog Number:

24562-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20089

GenBank Accession Number:

BC150311

GeneID (NCBI):

54662

UNIPROT ID:

Q9NVG8

Full Name:

TBC1 domain family, member 13

Calculated MW:

400 aa, 47 kDa

Observed MW:

47 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse heart tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24562-1-AP (TBC1D13 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.