

For Research Use Only

# GALNT13 Polyclonal antibody

Catalog Number: 24604-1-AP



## Basic Information

Catalog Number:

24604-1-AP

Size:

290 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16526

GenBank Accession Number:

BC101032

GeneID (NCBI):

114805

UNIPROT ID:

Q8IUC8

Full Name:

UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetyl galactosaminyltransferase 13 (GalNAc-T13)

Calculated MW:

561 aa, 65 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data