

For Research Use Only

GPR58/TAAR2 Polyclonal antibody

Catalog Number: 24609-1-AP



Basic Information

Catalog Number:

24609-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18939

GenBank Accession Number:

BC067461

GeneID (NCBI):

9287

UNIPROT ID:

Q9P1P5

Full Name:

trace amine associated receptor 2

Calculated MW:

351 aa, 40 kDa

Observed MW:

42 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : rat brain tissue, mouse brain tissue, mouse cerebellum tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat brain tissue were subjected to SDS PAGE followed by western blot with 24609-1-AP (GPR58/TAAR2 antibody at dilution of 1:600 incubated at room temperature for 1.5 hours.