

For Research Use Only

TMPRSS5 Polyclonal antibody

Catalog Number: 24616-1-AP



Basic Information

Catalog Number: 24616-1-AP	GenBank Accession Number: BC121803	Purification Method: Antigen affinity purification
Size: 400 µg/ml	GeneID (NCBI): 80975	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q9H3S3	
Isotype: IgG	Full Name: transmembrane protease, serine 5	
Immunogen Catalog Number: AG18493	Calculated MW: 457 aa, 50 kDa	
	Observed MW: 53 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : rat brain tissue,
Species Specificity: human, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat brain tissue were subjected to SDS PAGE followed by western blot with 24616-1-AP (TMPRSS5 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.