

For Research Use Only

B3GAT2 Polyclonal antibody

Catalog Number: 24631-1-AP



Basic Information

Catalog Number:

24631-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18402

GenBank Accession Number:

BC113995

GeneID (NCBI):

135152

UNIPROT ID:

Q9NPZ5

Full Name:

beta-1,3-glucuronyltransferase 2
(glucuronosyltransferase S)

Calculated MW:

323 aa, 37 kDa

Observed MW:

37 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse, rat

Positive Controls:

WB : mouse brain tissue, mouse cerebellum tissue,
mouse kidney tissue, mouse retina tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24631-1-AP (B3GAT2 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.