

For Research Use Only

MEP1A Polyclonal antibody, PBS Only

Catalog Number: 24640-1-PBS



Basic Information

Catalog Number:

24640-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18950

GenBank Accession Number:

BC131714

GeneID (NCBI):

4224

UNIPROT ID:

Q16819

Full Name:

meprin A, alpha (PABA peptide hydrolase)

Calculated MW:

746 aa, 84 kDa

Observed MW:

84 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24640-1-AP (MEP1A antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24640-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded colon cancer slide using 24640-1-AP (MEP1A antibody) at dilution of 1:300 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 24640-1-PBS in a different storage buffer formulation.