

For Research Use Only

GKN2 Polyclonal antibody

Catalog Number: 24651-1-AP



Basic Information

Catalog Number:

24651-1-AP

Size:

1200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19153

GenBank Accession Number:

BC146484

GeneID (NCBI):

200504

UNIPROT ID:

Q86XP6

Full Name:

gastrokine 2

Calculated MW:

184 aa, 20 kDa

Observed MW:

20 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse stomach tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse stomach tissue were subjected to SDS PAGE followed by western blot with 24651-1-AP (GKN2 Antibody) at dilution of 1:1000 incubated at 4 degree celsius over night.