

For Research Use Only

TAF13 Polyclonal antibody

Catalog Number: 24672-1-AP



Basic Information

Catalog Number:

24672-1-AP

Size:

1100 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18193

GenBank Accession Number:

BC121180

GeneID (NCBI):

6884

UNIPROT ID:

Q15543

Full Name:

TAF13 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 18kDa

Calculated MW:

124 aa, 14 kDa

Observed MW:

14-18 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse testis tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 24672-1-AP (TAF13 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.