

For Research Use Only

# TAF13 Polyclonal antibody

Catalog Number: 24672-1-AP



## Basic Information

Catalog Number:

24672-1-AP

Size:

1100 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18193

GenBank Accession Number:

BC121180

GeneID (NCBI):

6884

UNIPROT ID:

Q15543

Full Name:

TAF13 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 18kDa

Calculated MW:

124 aa, 14 kDa

Observed MW:

14-18 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse testis tissue,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 24672-1-AP (TAF13 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.