

For Research Use Only

ZNF284 Polyclonal antibody

Catalog Number: 24681-1-AP



Basic Information

Catalog Number:

24681-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20336

GenBank Accession Number:

BC156663

GeneID (NCBI):

342909

UNIPROT ID:

Q2VY69

Full Name:

zinc finger protein 284

Calculated MW:

593 aa, 69 kDa

Observed MW:

69 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse liver tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 24681-1-AP (ZNF284 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.