

For Research Use Only

ITSN2 Polyclonal antibody, PBS Only

Catalog Number: 24709-1-PBS



Basic Information

Catalog Number:

24709-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20130

GenBank Accession Number:

BC146779

GeneID (NCBI):

50618

UNIPROT ID:

Q9NZM3

Full Name:

intersectin 2

Calculated MW:

1697 aa, 193 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data