

For Research Use Only

ZFP2 Polyclonal antibody

Catalog Number: 24720-1-AP



Basic Information

Catalog Number:

24720-1-AP

Size:

5500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20087

GenBank Accession Number:

BC142989

GeneID (NCBI):

80108

UNIPROT ID:

Q6ZN57

Full Name:

zinc finger protein 2 homolog (mouse)

Calculated MW:

461 aa, 53 kDa

Observed MW:

53-55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HepG2 cells, Neuro-2a cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 24720-1-AP (ZFP2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.