

For Research Use Only

# MRPL47 Polyclonal antibody

Catalog Number: 24728-1-AP



## Basic Information

Catalog Number:

24728-1-AP

Size:

1000  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11744

GenBank Accession Number:

BC032522

GeneID (NCBI):

57129

UNIPROT ID:

Q9HD33

Full Name:

mitochondrial ribosomal protein L47

Calculated MW:

230 aa, 27 kDa

Observed MW:

29 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse testis tissue,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 24728-1-AP (MRPL47 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.