

For Research Use Only

GRIN3A Polyclonal antibody

Catalog Number: 24774-1-AP



Basic Information

Catalog Number:

24774-1-AP

Size:

340 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20416

GenBank Accession Number:

BC132866

GeneID (NCBI):

116443

UNIPROT ID:

Q8TCU5

Full Name:

glutamate receptor, ionotropic, N-methyl-D-aspartate 3A

Calculated MW:

1115 aa, 126 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data