

For Research Use Only

# SLC36A1 Polyclonal antibody, PBS Only

Catalog Number: 24775-1-PBS

Featured Product



## Basic Information

Catalog Number:

24775-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20594

GenBank Accession Number:

BC136437

GeneID (NCBI):

206358

UNIPROT ID:

Q7Z2H8

Full Name:

solute carrier family 36  
(proton/amino acid symporter),  
member 1

Calculated MW:

476 aa, 53 kDa

Observed MW:

43-53 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

## Background Information

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human small intestine tissue slide using 24775-1-AP (SLC36A1 antibody) at dilution of 1:200 (under 0x lens). This data was developed using the same antibody clone with 24775-1-PBS in a different storage buffer formulation.



HeLa cells were subjected to SDS PAGE followed by western blot with 24775-1-AP (SLC36A1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24775-1-PBS in a different storage buffer formulation.