

For Research Use Only

TMEM35 Polyclonal antibody

Catalog Number: 24786-1-AP

2 Publications



Basic Information

Catalog Number:

24786-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20589

GenBank Accession Number:

BC078658

GeneID (NCBI):

59353

UNIPROT ID:

Q53FP2

Full Name:

transmembrane protein 35

Calculated MW:

167 aa, 18 kDa

Observed MW:

20 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse cerebellum tissue, mouse brain tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yogesh Hooda	39553992	bioRxiv	
Niyati A Borkar	39236288	Am J Respir Cell Mol Biol	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse cerebellum tissue were subjected to SDS PAGE followed by western blot with 24786-1-AP (TMEM35 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.