

For Research Use Only

RXFP3 Polyclonal antibody, PBS Only

Catalog Number: 24798-1-PBS



Basic Information

Catalog Number:

24798-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20453

GenBank Accession Number:

BC113440

GeneID (NCBI):

51289

UNIPROT ID:

Q9NSD7

Full Name:

relaxin/insulin-like family peptide
receptor 3

Calculated MW:

469 aa, 51 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24798-1-AP (RFXP3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24798-1-PBS in a different storage buffer formulation.