

For Research Use Only

PRLHR Polyclonal antibody

Catalog Number: 24800-1-AP



Basic Information

Catalog Number:

24800-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20098

GenBank Accession Number:

BC101490

GeneID (NCBI):

2834

UNIPROT ID:

P49683

Full Name:

prolactin releasing hormone receptor

Calculated MW:

370 aa, 41 kDa

Observed MW:

41 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : HepG2 cells, Jurkat cells, SH-SY5Y cells, rat brain tissue

Background Information

PRLHR (prolactin-releasing hormone receptor) is implicated in lactation, regulation of food intake, and pain-signal processing.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24800-1-AP (PRLHR antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.