

For Research Use Only

POTEH Polyclonal antibody

Catalog Number: 24837-1-AP



Basic Information

Catalog Number:

24837-1-AP

Size:

280 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20275

GenBank Accession Number:

BC148758

GeneID (NCBI):

23784

UNIPROT ID:

Q6S545

Full Name:

POTE ankyrin domain family,
member H

Calculated MW:

545 aa, 61 kDa

Observed MW:

61 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : DU 145 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



DU 145 cells were subjected to SDS PAGE followed by western blot with 24837-1-AP (POTEH Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.