

For Research Use Only

C20orf152 Polyclonal antibody

Catalog Number: 24870-1-AP



Basic Information

Catalog Number:

24870-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21668

GenBank Accession Number:

BC114945

GeneID (NCBI):

140894

UNIPROT ID:

Q96M20

Full Name:

chromosome 20 open reading frame
152

Calculated MW:

576 aa, 68 kDa

Observed MW:

67-70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse lung tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse lung tissue were subjected to SDS PAGE followed by western blot with 24870-1-AP (C20orf152 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.