

For Research Use Only

DTWD2 Polyclonal antibody

Catalog Number: 24873-1-AP

1 Publications



Basic Information

Catalog Number:

24873-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21665

GenBank Accession Number:

BC114955

GeneID (NCBI):

285605

UNIPROT ID:

Q8NBA8

Full Name:

DTW domain containing 2

Calculated MW:

298 aa, 33 kDa

Observed MW:

33 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

IP: COLO 320 cells,

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|------------|-----------|---------|-------------|
| Yixuan Xie | 39173631 | Cell | WB |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-DTWD2 (IP:24873-1-AP, 3ug;
Detection:24873-1-AP 1:300) with COLO 320 cells
lysate 2000ug.