

For Research Use Only

URB2 Polyclonal antibody

Catalog Number: 24881-1-AP



Basic Information

Catalog Number:

24881-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21667

GenBank Accession Number:

BC114559

GeneID (NCBI):

9816

UNIPROT ID:

Q14146

Full Name:

URB2 ribosome biogenesis 2 homolog
(S. cerevisiae)

Calculated MW:

1524 aa, 171 kDa

Observed MW:

~150 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 24881-1-AP (URB2 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.