

For Research Use Only

TMEM190 Polyclonal antibody

Catalog Number: 24904-1-AP



Basic Information

Catalog Number: 24904-1-AP	GenBank Accession Number: BC128189	Purification Method: Antigen affinity purification
Size: 350 µg/ml	GeneID (NCBI): 147744	Recommended Dilutions: IHC 1:20-1:200
Source: Rabbit	UNIPROT ID: Q8WZ59	
Isotype: IgG	Full Name: transmembrane protein 190	
Immunogen Catalog Number: AG19131	Calculated MW: 177 aa, 19 kDa	

Applications

Tested Applications:
IHC, ELISA

Species Specificity:
human

Positive Controls:

IHC: human colon cancer tissue, human liver cancer tissue

Note-IHC: suggested antigen retrieval with **TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 24904-1-AP (TMEM190 Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 24904-1-AP (TMEM190 Antibody) at dilution of 1:50 (under 40x lens).