

For Research Use Only

IGF2 Polyclonal antibody

Catalog Number: 24911-1-AP

2 Publications



Basic Information

Catalog Number:

24911-1-AP

Size:

280 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21397

GenBank Accession Number:

BC000531

GeneID (NCBI):

3481

UNIPROT ID:

P01344

Full Name:

insulin-like growth factor 2
(somatomedin A)

Calculated MW:

236 aa, 26 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jinglin Zhang	36373776	J Pathol	IHC
Chao Zhang	38975094	Heliyon	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

