

For Research Use Only

ZNF114 Polyclonal antibody

Catalog Number: 24974-1-AP



Basic Information

Catalog Number:

24974-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21810

GenBank Accession Number:

BC125069

GeneID (NCBI):

163071

UNIPROT ID:

Q8NC26

Full Name:

zinc finger protein 114

Calculated MW:

417 aa, 48 kDa

Observed MW:

40-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : COLO 320 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



COLO 320 cells were subjected to SDS PAGE followed by western blot with 24974-1-AP (ZNF114 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.