

For Research Use Only

SPATA1 Polyclonal antibody

Catalog Number: 24980-1-AP



Basic Information

Catalog Number:

24980-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21784

GenBank Accession Number:

BC127634

GeneID (NCBI):

64173

UNIPROT ID:

Q5VX52

Full Name:

spermatogenesis associated 1

Calculated MW:

437 aa, 50 kDa

Observed MW:

40-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : human testis tissue, mouse testis tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 24980-1-AP (SPATA1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.