

For Research Use Only

ZNF586 Polyclonal antibody

Catalog Number: 25023-1-AP



Basic Information

Catalog Number:

25023-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21732

GenBank Accession Number:

BC124558

GeneID (NCBI):

54807

UNIPROT ID:

Q9NXT0

Full Name:

zinc finger protein 586

Calculated MW:

402 aa, 46 kDa

Observed MW:

46 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, COLO 320 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 25023-1-AP (ZNF586 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.