

For Research Use Only

CCDC40 Polyclonal antibody

Catalog Number: 25049-1-AP



Basic Information

Catalog Number:

25049-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20034

GenBank Accession Number:

BC035251

GeneID (NCBI):

55036

UNIPROT ID:

Q4G0X9

Full Name:

coiled-coil domain containing 40

Calculated MW:

1142 aa, 130 kDa

Observed MW:

90-100 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : fetal human brain tissue, A549 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



fetal human brain tissue were subjected to SDS PAGE followed by western blot with 25049-1-AP (CCDC40 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.