

For Research Use Only

# ARRDC3 Polyclonal antibody, PBS Only

Catalog Number: 25060-1-PBS

Featured Product



## Basic Information

Catalog Number:

25060-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18717

GenBank Accession Number:

BC053619

GeneID (NCBI):

57561

UNIPROT ID:

Q96B67

Full Name:

arrestin domain containing 3

Calculated MW:

414 aa, 46 kDa

Observed MW:

46 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data