

For Research Use Only

NYAP2 Polyclonal antibody

Catalog Number: 25108-1-AP



Basic Information

Catalog Number:

25108-1-AP

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19004

GenBank Accession Number:

BC131727

GeneID (NCBI):

57624

UNIPROT ID:

Q9P242

Full Name:

KIAA1486 protein

Calculated MW:

653 aa, 71 kDa

Observed MW:

80 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse brain tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 25108-1-AP (NYAP2) antibody at dilution of 1:1000 incubated at room temperature for 1.5 hours.