

For Research Use Only

# ARSI Polyclonal antibody, PBS Only

Catalog Number: 25112-1-PBS



## Basic Information

Catalog Number:

25112-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18884

GenBank Accession Number:

BC129995

GeneID (NCBI):

340075

UNIPROT ID:

Q5FYB1

Full Name:

arylsulfatase family, member I

Calculated MW:

569 aa, 64 kDa

Observed MW:

64 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 25112-1-AP (ARSI antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25112-1-PBS in a different storage buffer formulation.



IP result of anti-ARSI (IP:25112-1-AP, 4ug; Detection:25112-1-AP 1:500) with A549 cells lysate 800 ug. This data was developed using the same antibody clone with 25112-1-PBS in a different storage buffer formulation.