

FIBCD1 Polyclonal antibody

Catalog Number: 25125-1-AP

Basic Information

Catalog Number:

25125-1-AP

Size:

1400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18355

GenBank Accession Number:

BC032953

GeneID (NCBI):

84929

UNIPROT ID:

Q8N539

Full Name:

fibrinogen C domain containing 1

Calculated MW:

461 aa, 51 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:20-1:200

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC : human colon tissue,

Background Information

Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:


PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 25125-1-AP (FIBCD1 Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 25125-1-AP (FIBCD1 Antibody) at dilution of 1:50 (under 40x lens).