

For Research Use Only

# acetylated Tubulin (Lys40) Polyclonal antibody



Catalog Number: 25126-1-AP

2 Publications

## Basic Information

Catalog Number:

25126-1-AP

Size:

190  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_006009

GeneID (NCBI):

7846

UNIPROT ID:

Q71U36

Full Name:

tubulin, alpha 1a

Calculated MW:

52 kDa

Observed MW:

50-55 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

human

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Hidetoshi Tanioka	35861932	Jpn J Ophthalmol	IF
Luojia Liu	37949303	Biochim Biophys Acta Mol Cell Res	WB,IF

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data