

# C12orf54 Polyclonal antibody

Catalog Number: 25177-1-AP

## Basic Information

**Catalog Number:**

25177-1-AP

**Size:**

1000 µg/ml

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG15026

**GenBank Accession Number:**

BC031670

**GeneID (NCBI):**

121273

**UNIPROT ID:**

Q6X4T0

**Full Name:**

chromosome 12 open reading frame

54

**Calculated MW:**

127 aa, 14 kDa

**Purification Method:**

Antigen affinity purification

**Recommended Dilutions:**

IHC 1:20-1:200

## Applications

**Tested Applications:**

IHC, ELISA

**Species Specificity:**

human

**Note-IHC: suggested antigen retrieval with *TE buffer pH 9.0*; (\*) Alternatively, antigen retrieval may be performed with *citrate buffer pH 6.0***

**Positive Controls:**

IHC : human testis tissue, human brain tissue

## Background Information

### Storage

**Storage:**


Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**


PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 25177-1-AP (C12orf54 Antibody) at dilution of 1:50 (under 40x lens).



Immunohistochemical analysis of paraffin-embedded human brain tissue slide using 25177-1-AP (C12orf54 Antibody) at dilution of 1:50 (under 40x lens).