

For Research Use Only

JIP3 Polyclonal antibody

Catalog Number: 25212-1-AP

Featured Product

3 Publications



Basic Information

Catalog Number:

25212-1-AP

Concentration:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19561

GenBank Accession Number:

BC137123

GeneID (NCBI):

23162

UNIPROT ID:

Q9UPT6

Full Name:

mitogen-activated protein kinase 8
interacting protein 3

Calculated MW:

1336 aa, 147 kDa

Observed MW:

147 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, rat, mouse

Cited Species:

human, mouse

Positive Controls:

WB : rat brain tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Peiyuan Zhang	32427938	Nat Commun	WB
Hailin Liu	35693438	Saudi Pharm J	WB
Wei Zhang	40111412	JCI Insight	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat brain tissue were subjected to SDS PAGE followed by western blot with 25212-1-AP (JIP3 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.