

For Research Use Only

# LAMB4 Polyclonal antibody

Catalog Number: 25217-1-AP



## Basic Information

<b>Catalog Number:</b> 25217-1-AP	<b>GenBank Accession Number:</b> BC140804	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 900 µg/ml	<b>GeneID (NCBI):</b> 22798	<b>Recommended Dilutions:</b> WB 1:500-1:1000
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> A4D0S4	
<b>Isotype:</b> IgG	<b>Full Name:</b> laminin, beta 4	
<b>Immunogen Catalog Number:</b> AG19401	<b>Calculated MW:</b> 1761 aa, 194 kDa	
	<b>Observed MW:</b> 193 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : mouse skin tissue,
<b>Species Specificity:</b> human, rat, mouse	

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse skin tissue were subjected to SDS PAGE followed by western blot with 25217-1-AP (LAMB4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.