

For Research Use Only

# C15orf44 Polyclonal antibody

Catalog Number: 25227-1-AP **1 Publications**



## Basic Information

<b>Catalog Number:</b> 25227-1-AP	<b>GenBank Accession Number:</b> BC007991	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 260 µg/ml	<b>GeneID (NCBI):</b> 81556	<b>Recommended Dilutions:</b> WB 1:500-1:1000
<b>Source:</b> Rabbit	<b>Full Name:</b> chromosome 15 open reading frame 44	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 518 aa, 57 kDa	
<b>Immunogen Catalog Number:</b> AG19785	<b>Observed MW:</b> 57 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : HEK-293 cells, HeLa cells
<b>Cited Applications:</b> WB	
<b>Species Specificity:</b> human	
<b>Cited Species:</b> human	

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Sarah R Offley	36920904	Cell Rep	WB


## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25227-1-AP (C15orf44 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.