

For Research Use Only

C15orf44 Polyclonal antibody

Catalog Number: 25227-1-AP

1 Publications



Basic Information

Catalog Number:

25227-1-AP

Size:

260 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19785

GenBank Accession Number:

BC007991

GeneID (NCBI):

81556

Full Name:

chromosome 15 open reading frame

44

Calculated MW:

518 aa, 57 kDa

Observed MW:

57 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : HEK-293 cells, Hela cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Sarah R Offley	36920904	Cell Rep	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25227-1-AP (C15orf44 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.