

For Research Use Only

COMMD8 Polyclonal antibody

Catalog Number: 25237-1-AP **1 Publications**



Basic Information

Catalog Number: 25237-1-AP	GenBank Accession Number: BC008371	Purification Method: Antigen affinity purification
Size: 750 µg/ml	GeneID (NCBI): 54951	Recommended Dilutions: WB 1:500-1:1000 IF 1:20-1:200
Source: Rabbit	UNIPROT ID: Q9NX08	
Isotype: IgG	Full Name: COMM domain containing 8	
Immunogen Catalog Number: AG18712	Calculated MW: 183 aa, 21 kDa	
	Observed MW: 21 kDa	

Applications

Tested Applications: IF/ICC, WB, ELISA	Positive Controls: WB : HEK-293 cells, IF : HEK-293 cells,
Cited Applications: WB	
Species Specificity: human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Michael D Healy	37172566	Cell	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 25237-1-AP (COMMD8 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of HEK-293 cells using 25237-1-AP (COMMD8 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).