

For Research Use Only

FAM98C Polyclonal antibody

Catalog Number: 25274-1-AP



Basic Information

Catalog Number:

25274-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15042

GenBank Accession Number:

BC034383

GeneID (NCBI):

147965

UNIPROT ID:

Q17RN3

Full Name:

family with sequence similarity 98, member C

Calculated MW:

349 aa, 37 kDa

Observed MW:

37 kDa, 28 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : A431 cells,

Background Information

FAM98C has two isoforms with MW 28-30 kDa and 37-40 kDa.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 25274-1-AP (FAM98C Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.