

For Research Use Only

ASB14 Polyclonal antibody, PBS Only

Catalog Number:25300-1-PBS



Basic Information

Catalog Number:

25300-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18147

GenBank Accession Number:

BC114953

GeneID (NCBI):

142686

UNIPROT ID:

A6NK59

Full Name:

ankyrin repeat and SOCS box-
containing 14

Calculated MW:

587 aa, 65 kDa

Observed MW:

65 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25300-1-AP (ASB14 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 25300-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 25300-1-AP (ASB14 Antibody) at dilution of 1:50 (under 40x lens). This data was developed using the same antibody clone with 25300-1-PBS in a different storage buffer formulation.